

SUPPLEMENTARY TABLE

Supplementary Table 1. Primers used for Real-Time PCR.

Gene	Primer sequence	Accession number	Tann (°C)
Col1A1	sense 5'-tgttcctactcagccgtctgt-3' antisense 5'-gagccctcgcttccgtactc-3'	NM_007742.3	60
Col1A2	sense 5'-tgcccatctggtaaagaag-3' antisense 5'-accttgccacctgaacac-3'	NM_007743.2	60
TGFβ1	sense 5'-ggagcccgaagcggacta-3' antisense 5'-cgaatgtctgacgtattgaagaaca-3'	NM_011577.1	60
TGFβ RI	sense 5'-tagctgaaattgacctaattcctcg-3' antisense 5'-tgcggttatggcagatatagacc-3'	NM_009370.2	60
TGFβ RII	sense 5'-gggattgccatagctgtcat-3' antisense 5'-tgatggcacaattgtcactg-3'	NM_009371.3	60
MMP-2	sense 5'-agatcttcttctcaaggaccggtt-3' antisense 5'-ggctggtcagtggttgggta-3'	NM_008610.2	60
CRABP-II	sense 5'-tgatgaggaagatcgctgtg-3' antisense 5'-ttcactctcccattcacc -3'	NM_007759.2	60
GAPDH	sense 5'-aacttggcattgtggaagg-3' antisense 5'-cacattggggtaggaacac-3'	NM_001289726.1	60