

## SUPPLEMENTARY TABLES

**Supplementary Table 1. Primer used for QPCR analysis.**

<b>Name</b>	<b>Forward primer (5'-&gt;3')</b>	<b>Reverse primer (5'-&gt;3')</b>
hsa_circ_0001320	TGAGGAGCACCCCATAGC	CATCATAGCCACTGACACGG
FOXP1	AAGTTTGTCCCATGAGGAGC	CTGGAGGATCTGCTGCATTT
IL-6	ACTCACCTCTTCAGAACGAATTG	CCATCTTTGGAAGGTTTCAGGTTG
TNF-a	CCTCTCTCTAATCAGCCCTCTG	GAGGACCTGGGAGTAGATGAG
Col3a1	TGCTGCTGGTACTCCTGGTCTG	ACCTGGACCGCCTGGTTCAC
Acta2	GTGTTGCCCTGAAGAGCAT	GCTGGGACATTGAAAGTCTCA
Col1a1	GAGGGCCAAGACGAAGACATC	CAGATCACGTCATCGCACAAAC
CD44	CTGCCGCTTTGCAGGTGTA	CATTGTGGGCAAGGTGCTATT
Ki67	ACGCTGGTTACTATCAAAAGG	CAGACCCATTTACTTGTGTTGGA
GAPDH	TCGACAGTCAGCCGCATCTTCTTT	ACCAAATCCGTTGACTCCGACCTT

**Supplementary Table 2. Primer used for RNA pulldown.**

<b>Name</b>	<b>Positive-probe-biotin</b>	<b>Negative-probe-biotin</b>
hsa_circ_0001320	ACACGGGAACCTTAGAAATG	CATTTCTAAAGTTCCCGTGT